

Documentation for calculating relaxation time from NIST Polymerization Stress Tensometer (SRI 6005)

If you are using this software please cite:

Sarkar, S., Baker, P. J., Chan, E., Lin-Gibson, S., & Chiang, M. (2017). Quantifying the sensitivity of network structure and properties from simultaneous measurements during photopolymerization. *Soft Matter*.

Directory Organization:

We recommend that the code directory **get_tau** and experimental data folder be organized in the following way.

get_tau	(contains python code for calculating relaxation time)
GlobalDataFolder	(contains all experimental data)
ExperimentFolder	(contains all trials of specific experiment)
<TrialName>.txt	(text file containing the information about the trial setup)
<TrialName>.csv	(Results of experiment in CSV format)

Example <TrialName>.txt file

Experimental Conditions

Measurements: Stress, Temperature, Degree of conversion

Beam type: Rigid steel

Sample Description: 70_30_50%

Diameter of sample (cm): 0.2500

Sample height (cm): 0.2000

C-Factor: 0.6250

Position of sample along the beam (cm): 20.0000

Compliance of the beam (cm/N): 1.2322E-3

Integration time (ms): 25.0000

Scans to average: 3.0000

Sample rate (s): 0.2000

Irradiation Mode: Standard

Start time of irradiation (s): 20.0000

Light Intensity (mW/cm²): 500.0000

Duration of irradiation (s): 40.0000

Duration of experiment: 10.00 minutes and 0.00 seconds

Additional comments on the experiment:

The NIST tensometer records the experimental condition in the above format. If the user's instrument records the information in a different way, it must be converted into the above one.

Example <TrialName>.csv

Time (s)	Voltage (V)	Disp. (cm)	Stress (MPa)	DC (%)	deltaT
0.142216	0.00008416	0.00000107	0.00010018	2.926763	-0.12565
0.342656	0.00009116	0.00000107	0.00010918	1.078467	0.016231
0.543096	0.00026116	0.00000207	0.00031118	0.381026	0.300331
0.743536	0.00071016	0.00000707	0.00084418	2.357989	0.959891
0.943976	0.00145116	0.00001407	0.00172418	1.337635	2.112882

Lastly, we need a provide some additional information regarding for the material, which are provided using <inputfile>.txt.

Inputs for folder location

Superfolder = GlobalDataFolder

Subfolder = ExperimentFolder

Testname = TrialName

Information about the composite

Initial viscosity (in Pa.s) = 1.0

High-frequency Youngs modulus (GPa) = 4.5

Maximum shrinkage strain = 0.10

Coefficient of thermal expansion (in 1/T) = 1E-4,0.5E-4

Molecular density (in moles/m³) = 2500

Numerical calculation details

Calculation time (s) = 250.0

Regularization parameter = 1E-2

Denoising window = 5

```
# Analysis outname  
Casename = TrialAnalysis
```

The coefficient of thermal expansion asked for are the limiting values for the resin and the cured polymer at a high DC.

Once, the input file and rest of the information has been organized. The code can be executed as,

```
python get_tau/get_tau.py -i <inputfile>.txt
```

The outputs of the calculation will be stored in a folder <Casename> under the directory ExperimentFolder.